

Human Semaphorin-4D/SEMA4D/CD100 protein (Recombinant) (C-His) (STJP001081)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10507
Gene Symbol	SEMA4D
UniProt ID	SEMA4D_HUMAN
Immunogen Region	Met22-Arg734
Immunogen Sequence	MAFAPIPRITWEHREVHLVQ FHEPDIYNYSALLSEDKDT LYIGAREAVFAVNALNISEK QHEVYWKVSEDKKAKCAEKG KSKQTECLNYIRVLQPLSAT SLYVCGTNAFQACDHLNLT SFKFLGKNEDGKGRCPFDPA HSYTSVMVDGELYSGTSYNF LGSEPIISRNSHSPRLRTEY AIPWLNPSFVFADVIRKSP DSPGEGDDRVYFFFEVSVE YEFVFRVLIPRIARVCKGD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

