

Human Follistatin-like protein 3/FSTL3 protein (Recombinant) (C-His) (STJP001086)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10272
Gene Symbol	FSTL3
UniProt ID	FSTL3_HUMAN
Immunogen Region	Met27-Val263
Immunogen Sequence	MGSGNPAPGGVCWLQQGQEA TCSLVLQTDVTRAECASGN IDTAWSNLTHPGNKINLLGF LGLVHCLPCKDSCDGVCEGP GKACRMLGGRPRCECAPDCS GLPARLQVCGSDGATYRDEC ELRAARCRGHPDLSVMYRGR CRKSCEHVVCPRPQSCVVDQ TGAHCVVCRAAPCPVPSSP GQELCGNNNVTYISSCHMRQ ATCFLGRSIGVRHAGSCAGT PEEPPGGESAEEEENFV

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

