

Human Tetranectin/CLEC3B protein (Recombinant) (C-His) (STJP001088)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7123
Gene Symbol	CLEC3B
UniProt ID	TETN_HUMAN
Immunogen Region	Glu22-Val202
Immunogen Sequence	EPPTQKPKKIVNAKKDVVNT KMFEELKSRLDTLAQEVALL KEQQALQTVCLKGTKVHMKC FLAFTQTKTFHEASEDCISR GGTLSTPQTGSENDALYEYL RQSVGNEAEIWLGLNDMAAE GTWVDMTGARIAYNWETEI TAQPDGGKTENCAVLGSAAN GKWFDKRCRDQLPYICQFGI V

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

