

Human CD14 protein (Recombinant) (C-His) (STJP001089)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	929
Gene Symbol	CD14
UniProt ID	CD14_HUMAN
Immunogen Region	Thr20-Cys352
Immunogen Sequence	TTPEPCELDDEDFRCVCNFS EPQPDWSEAFQCVSAVEVEI HAGGLNLEPFLKRVDADADP RQYADTVKALRVRLTVGAA QVPAQLLVGALRVLAYSRK ELTLEDLKITGTMPPPLEA TGLALSSLRLRNVSATGRS WLAELQQWLKPKLVLSIAQ AHSPAFSCEQVRAFPALTS DLSDNPGLGERGLMAALCPH KFPAIQNLALRNTGMETPTG VCAALAAAGVQPHSLDLSH

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

