

Human CDH5/VE-cadherin protein (Recombinant) (C-hFc&His) (STJP001090)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1003
Gene Symbol	CDH5
UniProt ID	CADH5_HUMAN
Immunogen Region	Asp48-Gln593
Immunogen Sequence	DWIWNQMHI DEEKNTSLPHH VGKIKSSVSRKNAKYLLKGE YVGKVFRVDAETGDVFAIER LDRENISEYHLTAVIVDKDT GENLETPSSFTIKVHDVNDN WPVFTHRLFNASVPESAVG TSVISVTAVDADDPTVGDHA SVMYQILKGKEYFAIDNSGR IITITKSLDREKQARYEIVV EARDAQGLRGDSGTATVLVT LQDINDNFPFFTQTKYTFVV PEDTRVGTSGSLFVEDPD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

