

Human Angiotensinogen/AGT protein (Recombinant) (C-His) (STJP001103)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	183
Gene Symbol	AGT
UniProt ID	ANGT_HUMAN
Immunogen Region	Asp 34-Ala 485
Immunogen Sequence	DRVYIHPFHLVIHNESTCEQ LAKANAGKPKDPTFIPAPIQ AKTSPVDEKALQDQLVLVAA KLDTEKDLRAAMVGMLANFL GFRIYGMHSELWGTVVHGATV LSPTAVFGTLASLYLGALDH TADRLQAILGVPWKDKNCTS RLDAAHKVLSALQAVQGLLVA QGRADSQAQLLLSTVVGFT APGLHLKQPFVQGLALYTPV VLPRSLDFTELDVAEEKIDR FMQAVTGWKTGCSLTGASV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

