

Human IL-13RA1/CD213a1 protein (Recombinant) (C-hFc&His) (STJP001106)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3597
Gene Symbol	IL13RA1
UniProt ID	I13R1_HUMAN
Immunogen Region	Met 1-Thr 343
Immunogen Sequence	MEWPARLCGLWALLLCAGGG GGGGGAAPTETQPPVTNLSV SVENLC TVIWTWNPPEGASS NCSLWYF SHFGDKQDKKIAP ETRRSIEVPLNERICLQVGS QCSTNESEKPSILVEKICSP PEGDPESA VTELCQIWHNLS YMKCSWLPGRNTSPDTNYTL YYWHR SLEKIHQCENIFREG QYFGCSFDLTKVKDSSFEQH SVQIMVKDNAGKIKPSFNIV PLTSRVKPDPPHIKNLSFH

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

