

Human MSR1/CD204 protein (Recombinant) (C-His) (STJP001113)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4481
Gene Symbol	MSR1
UniProt ID	MSRE_HUMAN
Immunogen Region	Lys77-Leu451
Immunogen Sequence	KWETKNCVSSTNANDITQS LTGKGNDSSEEMRFQEVFME HMSNMEKRIQHILDMEANLM DTEHFQNFSTMTDQRFNDIL LQLSTLFSSVQGHGNAIDEI SKSLISLNTTLLDLQLNIEN LNGKIQENTFKQQEEISKLE ERVYNVSAEIMAMKEEQVHL EQEIKGEVKVLNNITNDLRL KDWEHSQTLRNITLIQGPPG PPGKGDGRGPTGESGPRGFP GPIGPPGLKGDRAIGFPG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

