

Human FcRn protein (Active Recombinant) (C-His (FCGRT) &No tag (B2M)) (STJP001116)

GENERAL INFORMATION

Product Type Proteins
Host / Source HEK293 cells

PRODUCT PROPERTIES

Formulation Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region Ile21-Met119
Immunogen Sequence AESHLSLLYLHTAVSSPAPG TPAFWVSGWLGPQQYLSYNS LRGEAEPCGAWVWENQVSWY
WEKETDDLRIKEKLFLEAFK ALGGKGPYTLQGLLGCELGP DNTSVPTAKFALNGEEFMNF
DLKQGTWGGDWPEALAIQR WQQQDKAANKELTFLFSCP HRLREHLERGRGNLEWKEPP
SMRLKARPSSPGFSVLTCSA FSFYPELQLRFLRNGLAAG TQQGDFGPNSDGFSFHASS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

