

Human IFN-gamma protein (Active Recombinant) (No tag) (STJP001117)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3458
Gene Symbol	IFNG
UniProt ID	IFNG_HUMAN
Immunogen Region	Gln24-Gln166
Immunogen Sequence	QDPYVKEAENLKKYFNAGHS DVADNGTLFLGILKNWKEES DRKIMQSQIVSFYFKLFKNF KDDQSIQKSVETIKEDMNVK FFNSNKKKRDDFEKLTNYSV TDLNVQRKAIHELIQVMAEL SPAAGTGKRKRSQLFRGRR ASQ

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

