

Human IL-7 protein (Recombinant) (C-His) (STJP001120)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3574
Gene Symbol	IL7
UniProt ID	IL7_HUMAN
Immunogen Region	Asp26-His177
Immunogen Sequence	DCDIEGKDGKQYESVLMVSI DQLLDSMKEIGSNCLNNEFN FFKRHICDANKEGMFLFRAA RKLRQFLKMNSTGDFDLHLL KVSEGTILLNCTGQVKGRK PAALGAEAQPTKSLEENKSLK EQKLLNDLCLFLKRLLEIKT CWNKILMGTKHE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

