

Human FGF-2/bFGF protein (Recombinant) (No tag) (STJP001122)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM Tris, 150 mM NaCl, pH 7.5.
Molecular Weight	16.4 kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2247
Gene Symbol	FGF2
UniProt ID	FGF2_HUMAN
Immunogen Region	Pro143-Ser288
Immunogen Sequence	PALPEDGGSGAFPPGHFKDP KRLYCKNGGFFLRHPDGRV DGVREKSDPHIKLQLQAEER GVVSIKGVCANRYLAMKEDG RLLASKCVTDECFFFERLES NNYNTYRSRKYTSWYVALKR TGQYKLGSKTGPGQKAILFL PMSAKS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

