

Human BTLA protein (Active Recombinant) (C-hFc&His) (STJP001130)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	151888
Gene Symbol	BTLA
UniProt ID	BTLA_HUMAN
Immunogen Region	Lys31-Thr134
Immunogen Sequence	KESC DVQLYIKRQSEHSILA GDPFELECPVKYCANRPHVT WCKLNGTTCVKLEDRQTSWK EEKNISFFILHFEPVLPNDN GSYRCSANFQSNLIESHSTT LYVT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

