

Mouse Flt3 ligand/Flt3L protein (Recombinant) (C-hFc&His) (STJP001138)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	14256
Gene Symbol	Flt3lg
UniProt ID	FLT3L_MOUSE
Immunogen Region	Met1-Arg188
Immunogen Sequence	MTVLPAWSPNSSLLLLLLL LSPCLRGTPDCYFHSPISS NFKVKFRELDHLLKDYPVT VAVNLQDEKHCKALWSLFLA QRWIEQLKTVAGSKMQTLLE DVNTEIHFTSCTFQPLPEC LRFVQTNISHLLKDTCTQLL ALKPCIGKACQNFRCLEVQ CQPDSSSTLLPPRSPIALEAT ELPEPRPR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

