

Mouse Dkk-1 protein (Recombinant) (C-His) (STJP001139)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13380
Gene Symbol	Dkk1
UniProt ID	DKK1_MOUSE
Immunogen Region	Met1-His272
Immunogen Sequence	MMVVCAAAVRFLAVFTMMA LCSLPLL GASATLNSVLINS NAIKNLPPPLGGAGGQPGSA VSVAPGVLYEGGNKYQLDN YQPYPCAEEDEECGSDEYCSS PSRGAAGVGGVQICLACRKR RKRCMRHAMCCPGNYCKNGI CMPSDHSHFPRGEIEESIIE NLGNDHNAAAGDGYPRRTL TSKIYHTKGQEGSVCLRSSD CAAGLCCARHFWSKICKPVL KEGQVCTKHKRKGSHGLEI

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

