

Human CD6 protein (Recombinant) (C-hFc&His) (STJP001160)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	923
Gene Symbol	CD6
UniProt ID	CD6_HUMAN
Immunogen Region	His18-Glu398
Immunogen Sequence	HPSPAPPDQLNTSSAESELW EPGERLPVRLTNGSSSCSGT VEVRLASWEPACGALWDSR AAEAVCRALGCGGAEASQL APPTPELPPPPAAGNTSVAA NATLAGAPALLCSGAEWRLC EVVEHACRS DGRRARVTCAE NRALRLVDGGACAGRVEML EHGEWGSVCDDTDWLEDAHV VCRQLGCGWAVQALPGLHFT PGRGPIHRDQVNCSGAEAYL WDCPGLPGQHYCGHKEDAG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

