

Mouse B7-H1/PD-L1/CD274 protein (Recombinant) (C-hFc&His) (STJP001164)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	60533
Gene Symbol	Cd274
UniProt ID	PD1L1_MOUSE
Immunogen Region	Phe19-Thr238
Immunogen Sequence	FTITAPKDLVYVVEYGSNTM ECRFPVERELDLLALVYWE KEDEQVIQFVAGEEDLKPQH SNFRGRASLPKDQLLKGNAALQITDVKLQDAGVYCCIISY GGADYKRITLKVNPYRKIN QRISVDPATSEHELICQAEG YPEAEVIWTNSDHQPVSGKR SVTTSRTEGMLLNVTSSLRV NATANDVFYCTFWRSQPGQN HTAELIPELPAHPQNR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

