

## Human VEGF-A/VEGF165 protein (Recombinant) (N-His) (STJP001165)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">7422</a>
Gene Symbol	<a href="#">VEGFA</a>
UniProt ID	<a href="#">VEGFA_HUMAN</a>
Immunogen Region	Ala27-Arg191
Immunogen Sequence	APMAEGGGQNHHEVVKFMDV YQRSYCHPIETLVDFQEYP DEIEYIFKPSCVPLMRCGGC CNDEGLECVPTTEESNITMQI MRIKPHQGGQHIGEMSFLQHN KCECRPKKDRARQENPCGPG SERRKHLFVQDPQTKCCK NTDSRCKARQLELNERTCRC DKPRR

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

