

Mouse TIM-3/HAVCR2/CD366 protein (Recombinant) (C-hFc&His) (STJP001169)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	171285
Gene Symbol	Havcr2
UniProt ID	HAVR2_MOUSE
Immunogen Region	Arg20-Arg191
Immunogen Sequence	RSLENAYVFEVGKNAYLPCS YTLSTPGALVPMCWGKGFCEP WSQCTNELLRTDERNVITYQK SSRYQLKGDNLNKGDVSLIK NVTLDHGTCCRIQFPGLM NDKKLEKLDIAAKVTPAQ TAHGDSTTASPRTLTERNG SETQTLVTLHNNNGTKISTW ADEIKDSGETIR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

