

Mouse ErbB-2/HER2/CD340 protein (Recombinant) (C-hFc&His) (STJP001182)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13866
Gene Symbol	ErbB2
UniProt ID	ERBB2_MOUSE
Immunogen Region	Thr23-Thr653
Immunogen Sequence	TQVCTGTDMLRLPASPETH LDMRLHLYQGCQVVGQNLLEL TYLPANASLSFLQDIQEVQG YMLIAHNRVKHVPLQRLRIV RGTQLFEDKYALAVLDNRDP LDNVTTAAPGRTPPEGLRELG LRSLTEILKGGVLIRGNPQL CYQDMVLWKDVLKNNQLAP VMDMTNRSRACPPCAPTCKD NHCWGESPEDCQILTGTICT SGCARCKGRLPTDCHEQCA AGCTGPKHSDCLACLHFNH

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

