

Human Annexin A5/ANXA5 protein (Recombinant) (C-His) (STJP001194)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	308
Gene Symbol	ANXA5
UniProt ID	ANXA5_HUMAN
Immunogen Region	Met1-Asp320
Immunogen Sequence	MAQVLRGTVTDFPGFDERAD AETLRKAMKGLGTDEESILT LLTSRSNAQRQEISAFAKTL FGRDLLDDLKSELTGKFEKL IVALMKPSRLYDAYELKHAL KGAGTNEKVLTEIIASRTPE ELRAIKQVYEEYGSLEDD VVGDTSGYYQRLVLLQAN RDPDAGIDEAQVEQDAQALF QAGELKWGTDEEKFITIFGT RSVSHLRKVFDKYMISGFQ IETIDRETSGNLEQLLLA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

