

Human Endothelial protein C receptor/CD201 protein (Recombinant) (C-His) (STJP001195)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10544
Gene Symbol	PROCR
UniProt ID	EPCR_HUMAN
Immunogen Region	Ser18-Ser210
Immunogen Sequence	SQDASDGLQLHMLQISYFR DPYHVWYQGNASLGGHLTHV LEGPDTNTTIIQLQLQEPE SWARTQSGLSYLLQFHGLV RLVHQERTLAFPLTIRCFLG CELPPEGSRAHVFFEAVNG SSFVSFRPERALWQADTQVT SGVVTFTLQQLNAYNRTRYE LREFLEDTCVQYVQKHISAE NTKGSQTSRSYTS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

