

Human LILRB4/ILT-3/CD85k protein (Recombinant) (C-His) (STJP001196)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11006
Gene Symbol	LILRB4
UniProt ID	LILRB4_HUMAN
Immunogen Region	Gln22-Glu259
Immunogen Sequence	QAGPLPKPTLWAEPPGVSISW GNSVTIWCQGTLEAREYRLD KEESPAPWDRQNPLEPKNKA RFSIPSMTEYAGRYRCYYR SPVGWSQPSDPLELVMGTAY SKPTLSALPSPLVTSGKSVT LLCQSRSPMDTFLLIKERAA HPLLHLRSEHGAQQHQAEFP MSPVTSVHGGTYRCFSSHGF SHYLLSHSPDPLELIVSGSL EGPRPSPTRSVSTAAGPEDQ PLMPTGSPVPHSGLRRHWE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

