

Human LYPD3 protein (Recombinant) (C-His) (STJP001201)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27076
Gene Symbol	LYPD3
UniProt ID	LYPD3_HUMAN
Immunogen Region	Leu31-His286
Immunogen Sequence	LECYSCVQKADDGCSPNKMK TVKCAPGVDVCTEAVGAVET IHGQFSLAVRGCGLPGKN DRGLDLHGLLAFIQLQQCAQ DRCNAKLNLTSRALDPAGNE SAYPPNGVECYSCVGLSREA CQGTSPVVSVCYNASDHVYK GCFDGNVTLTAANVTVSLPV RGCYQDEFCTR DGVTGPGFT LSGSCCQGSRCNSDLRNKTY FSPRIPLVRLPPPEPTTVA STTSVTTSTSAVPRPTSTT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

