

Human NovH/CCN3 (N97K) protein (Recombinant) (C-His) (STJP001202)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4856
Gene Symbol	CCN3
UniProt ID	CCN3_HUMAN
Immunogen Region	Thr32-Met357
Immunogen Sequence	TQRCPPQCPGRCPATPPTCA PGVRAVLGDCSCCLVCARQR GESCSLDLEPCDESSGLYCDR SADPSKQTGICTAVEGDNCV FDGVIYRSGEKFQPSCKFQC TCRDGGIGCVPRCQLDVLLP EPNCPAPRKVEVPGECECKW ICGPDEEDSLGGLTLAAYRP EATLGVEVSDSSVNCIEQTT EWTACSKSCGMGFSTRVTNR NRQCEMLKQTRLCMVRPCEQ EPEQPTDKKGKKCLRTKKS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

