

Human MERTK protein (Active Recombinant) (C-hFc&His) (STJP001203)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10461
Gene Symbol	MERTK
UniProt ID	MERTK_HUMAN
Immunogen Region	Arg26-Ala499
Immunogen Sequence	REEAKPYPLFPGPFPGSLQT DHTPLLSLPHASGYQPALMF SPTQPGRPHTGNVAIPQVTS VESKPLPPLAFKHTVGHIL SEHKGVKFNCSISVPNIYQD TTISWWKDGKELLGAHHAIT QFYPDDEVTAIASFISITSV QRSDNGSYICKMKINNEEIV SDPIYIEVQGLPHFTKQPES MNVTRNTAFNLTCQAVGPPE PVNIFVWQNSSRVNEQPEKS PSLVTVPLGLETEMAVFSCEA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

