

Human FLT-1/VEGFR-1 protein (Recombinant) (C-His) (STJP001205)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2321
Gene Symbol	FLT1
UniProt ID	VGFR1_HUMAN
Immunogen Region	Met1-Asn756
Immunogen Sequence	MVSYWDTGVLLCALLSCLLL TGSSSGSKLKDPPELSLKGTO HIMQAGQTLHLQCRGEAAHK WSLPEMVSKESERLSITKSA CGRNGKQFCSTLTLNTAQAN HTGFYSCKYLAVPTSKKKET ESAIYIFISDTGRPFVEMYS EIPEIIHMTGRELVIPCRV TSPNITVTLKKFPLDTLIPD GKRIIWDSRKGFIISNATYK EIGLLTCEATVNGHLYKTNV LTHRQTNTIIVQISTPRP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

