

Human Siglec-9/CD329 protein (Recombinant) (C-hFc&His) (STJP001208)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27180
Gene Symbol	SIGLEC9
UniProt ID	SIGL9_HUMAN
Immunogen Region	Gln18-Gly348
Immunogen Sequence	QTSKLLTMQSSVTVQEGLCV HVPCSFSSYP SHGWYIPGPV V HGYWVFREGANTDQDAPVATN NPARAVWEETRDRFHLGDP HTKNCTLSIRDARRSDAGRY FFRMEKGSIKWNYKHHRLSV NVTALTRHPNIIIPGTLESG CPQNLTCVSPWACEQGT PPM ISWIGTSVSP LDPSTTRSSV LTLIPQPQDHGTS LTCQVTF PGASVTTNKT VHLNVSYP PQ NLTMTVFQGDGT VSTVLGN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

