

Human LRRC4 protein (Recombinant) (C-His) (STJP001209)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64101
Gene Symbol	LRRC4
UniProt ID	LRRC4_HUMAN
Immunogen Region	Ala39-Lys527
Immunogen Sequence	ASAGPQNCPSVCSCSNQFSK VVCTRRGLSEVPQGIPSNTR YLNLMEINNIQMIQADTFRHL HHLEVLQLGRNSIRQIEVGA FNGLASLNTLELFDNWLTVI PSGAFEYLSKLRRELWLRNNP IESIPSYAFNRVPSLMRLDL GELKKLEYISEGAFEGFLNL KYLNLGMCNIKDMPNLTPLV GLEELEMSGNHFPPIRPGSF HGLSSLKKLWVMNSQVSLIE RNAFDGLASLVELNLAHNN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

