

Human Layilin/LAYN protein (Recombinant) (C-hFc&His) (STJP001210)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	143903
Gene Symbol	LAYN
UniProt ID	LAYN_HUMAN
Immunogen Region	Met1-Glu220
Immunogen Sequence	MRPGTALQAVLLAVLLVGLR AATGRLLSGQPVCRGGTQRP CYKVIYFHDTSRRLNFEEAK EACRRDGGQLVSIESTEDEQK LIEKFIENTLLPSDGDWIGL RRREEKQSNSTACQDLYAWT DGSISQFRNWYVDEPSCGSE VCVVMYHQPSAPAGIGGPYM FQWNDDRCNMKNNFICKYSD EKPAVPSREAEGEETELTTP VLPEETQEEDAKKTFKESRE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

