

Human Siglec-2/CD22 protein (Recombinant) (C-hFc&His) (STJP001214)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	933
Gene Symbol	CD22
UniProt ID	CD22_HUMAN
Immunogen Region	Asp20-Arg687
Immunogen Sequence	DSSKWFVPEHPETLYAWEGAC VWIPCTYRALDGDLESFILF HNPEYNKNTSKFDGTRLYES TKDGKVPSEQKRVQFLGDKN KNCTLSIHPVHLNDSGQLGL RMESKTEKWMERIHNLNVSER PFPPIQLPPEIQESQEVTL TCLLNFSYGYPIQLQWLLE GVPMRQAAVTSTSLTIKSVF TRSEKLFSPQWSHHGKIVTC QLQDADGKFLSNDTVQLNVK HTPKLEIKVTPSDAIVREG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

