

Human CADM1/IGSF4A protein (Recombinant) (C-His) (STJP001215)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23705
Gene Symbol	CADM1
UniProt ID	CADM1_HUMAN
Immunogen Region	Gln45-His374
Immunogen Sequence	QNLFTKDVTVIEGEVATISC QVNKSDDSVIQLLNPNRQTI YFRDFRPLKDSRFQLNFSS SELKVSLTNVSISDEGRYFC QLYTDPPQESYTTITVLVPP RNLMDIQKDTAVEGEEIEV NCTAMASKPATTIRWFKGNT ELKGKSEVEEWSMDMYTVTSQ LMLKVHKEDDGVVICQVEH PAVTGNLQTQRYLEVQYKPQ VHIQMTYPLQGLTREGDALE LTCEAIGKPQPMVTWVRV

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

