

Human ADAM8/CD156a protein (Recombinant) (C-His) (STJP001216)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M filtered solution of 50mM Tris, 10mM CaCl ₂ , 150mM NaCl, pH 7.5.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	101
Gene Symbol	ADAM8
UniProt ID	ADAM8_HUMAN
Immunogen Region	Met1-Pro497
Immunogen Sequence	MRGLGLWLLGAMMLPAIAPS RPWALMEQYEVLPWRLPGP RVRRALPSHLGLHPERVSYV LGATGHNFTLHLRKNRDLLG SGYTETYTAANGSEVTEQPR GQDHCYQGHVEGYPDASAAS LSTCAGLRGFFQVGSDDLHI EPLDEGGEGGRHAVYQAEHL LQTAGTCGVSDDSLGSLLGP RTAAVFRPRPGDSLPSRETR YVELYVVVDNAEFQMLGSEA AVRHRVLEVNVHVDKLYQK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

