

Human Frizzled-5/FZD5 protein (Recombinant) (C-His) (STJP001232)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7855
Gene Symbol	FZD5
UniProt ID	FZD5_HUMAN
Immunogen Region	Met1-Pro167
Immunogen Sequence	MARPDPSAPPSLLLLLAQL VGRAAAASKAPVCQEITVPM CRGIGYNLTHMPNQFNHDTQ DEAGLEVHQFWPLVEIQCSPLDLRFFLCSMYTPICLPDYHK PLPPCRSV CERAKAGCSPLM RQYGF AWP PERMSCDRLPVLGRDAEVL CMDYNRSEATTAPP RPFPAKP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

