

Human Noggin/NOG protein (Recombinant) (C-His) (STJP001239)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES







Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9241
Gene Symbol	NOG
UniProt ID	NOGG_HUMAN
Immunogen Region	Gln28-Cys232
Immunogen Sequence	QHYLHIRPAPSDNLPVLDLI EHPDPIFDPKKDLNETLLR SLLGGHYDPGFMATSPPEDR PGGGGGAAGGAEDLAELDQL LRQRPSGAMPSEIKGLEFSE GLAQQGKQRLSKLRRKLQM WLWSQTFCPVLYAWNDLGSR FWPRYVKVGSCFSKRSCSVP EGMVCKPSKSVHLTVLRWRC QRRGGQRCGWIPQIPIISE CKCSC

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



170kDa —
 130kDa —
 100kDa —
 70kDa —
 55kDa —
 40kDa —
 35kDa —
 25kDa —
 15kDa —
 10kDa —