

Human Fc-gamma RI/CD64 protein (Recombinant) (C-His) (STJP001241)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2209
Gene Symbol	FCGR1A
UniProt ID	FCGR1_HUMAN
Immunogen Region	Met1-Pro288
Immunogen Sequence	MWFLTLLLWVPVDGQVDTT KAVITLQPPWVSVFQEETVT LHCEVLHLPGSSSTQWFLNG TATQTSTPSYRITSASVNDG GEYRCQRGLSGRSDPIQLEI HRGWLLQLVSSRVFTEGEPL ALRCHAWKDKLVYNVLYYRN GKAFKFFHWNSNLTKTNI SHNGTYHCSGMGKHRYTSAG ISVTVKELFPAPVLNASVTS PLLEGNLVTLSCKETKLLLR PGLQLYFSFYMGSKTLRGR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

