

Human CD24 protein (Active Recombinant) (C-hFc) (STJP001243)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100133941
Gene Symbol	CD24
UniProt ID	CD24_HUMAN
Immunogen Region	Met1-Gly59
Immunogen Sequence	MGRAMVARLGLGLLLALLL PTQIYSSETTTGTSSNSSQS TSNSGLAPNPTNATTKAAG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

