

## Human FGF-19 protein (Recombinant) (C-His) (STJP001244)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">9965</a>
Gene Symbol	<a href="#">FGF19</a>
UniProt ID	<a href="#">FGF19_HUMAN</a>
Immunogen Region	Phe27-Lys216
Immunogen Sequence	FSDAGPHVHYGWDPPIRLRH LYTSGPHGLSSCFLRIRADG VVDCARGQSAHSLLEIKAVA LRTVAIKGVHSVRYLCLMGAD GKMQLLQYSEEDCAFEED RPDGYNVYRSEKHRLPVSL SAKQRQLYKNRGLPLSHFL PMLPMVPEEPEDLRGHLESD MFSSPLETDSMDPFGLVTGL EAVRSPSFEK

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

