

Mouse B7-1/CD80 protein (Recombinant) (C-His) (STJP001245)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12519
Gene Symbol	Cd80
UniProt ID	CD80_MOUSE
Immunogen Region	Val38-Lys245
Immunogen Sequence	VDEQLSKSVKDKVLLPCRYN SPHEDESEDRIYWQKHDKVV LSVIAGKLVWPEYKNRTLY DNTTYSLIILGLVLSDRGTY SCVVQKKERGTYEVKHLALV KLSIKADFSTPNITESGNPS ADTKRITCFASGGFPKPRFS WLENGRELPGINTTISQDPE SELYTISSQLDFNTTRNHTI KCLIKYGAHVSEDFTWKEP PEDPPDSK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

