

## Human IL-15 protein (Active Recombinant) (C-His) (STJP001250)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.                                 |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                              |
|--------------------|------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">3600</a>                                                                                                         |
| Gene Symbol        | <a href="#">IL15</a>                                                                                                         |
| UniProt ID         | <a href="#">IL15_HUMAN</a>                                                                                                   |
| Immunogen Region   | Asn49-Ser162                                                                                                                 |
| Immunogen Sequence | NWVNVISDLKKIEDLIQSMH IDATLYTESDVHPSCCKVTAM KCFLLELQVISLESQDASIH<br>DTVENLIILANNSLSSNGNV TESGCKECEEELEEKNIKEFL QSFVHIVQMFINTS |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

