

Human GTPase KRas/KRAS (G12V) protein (Recombinant) (N-His) (STJP001252)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3845
Gene Symbol	KRAS
UniProt ID	P04461
Immunogen Region	Thr2-Cys186
Immunogen Sequence	TEYKLVVVGAVGVGKSALTIQLIQNHFVDEYDPTIEDSYR KQVVIDGETCLLDILDITAGQ EEYSAMRDQYMRGTGEGFLCV FAINNTKSFEDIHYREQIK RVKDSQEDVPMVLVGNKCDLP SRTVDTKQAQDLARSYGIPF IETSAKTRQRVEDAFYTLVR EIRQYRLKKISKEEKTTPGCV KIKKC

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

