

Human BDNF protein (Active Recombinant) (C-His) (STJP001258)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	627
Gene Symbol	BDNF
UniProt ID	BDNF_HUMAN
Immunogen Region	His129-Arg247
Immunogen Sequence	HSDPARRGELSVCDISEWV TAADKKTAVDMSGGTVTVLE KVPVSKGQLKQFYFETKCNP MGYTKEGCRGIDKRHWNSQC RTTQSYVRALTMDSKKRIGW RFIRIDTSCVCTLTIKRGR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

