

Human IL-15 protein (Active Recombinant) (C-hFc&His) (STJP001264)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3600
Gene Symbol	IL15
UniProt ID	IL15_HUMAN
Immunogen Region	Asn49-Ser162
Immunogen Sequence	NWVNVISDLKKIEDLIQSMH IDATLYTESDVHPSCKVTAM KCFLLLELQVISLES GDASIH DTVENLILANNLSLSSNGNV TEGSCKECEELEEKNIKEFL QSFVHIVQMFINTS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

