

SARS-CoV-2 Spike S1 protein (Active Recombinant) (C-His) (STJP001268)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4. or Supplied as a 0.22 Mu m filtered solution in PBS, pH 7.4.
Molecular Weight	76kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Val11-Arg682
Immunogen Sequence	VSSQCVNLTTRTQLPPAYTN SFTRGVYYPDKVFRSSVLHS TQDLFLPFFSNVTWFHAIHV SGTNGTKRFDNPVLPFNDGV YFASTKSNIRGWIFGTTL DSKTQSLIVNNATNVVIKV CEFOFCNDPFLGVYHKNNK SWMESEFRVYSSANNCTFEY VSQPFLLMDLEGKQGNFKNLR EFVFKNIDGYFKIYSKHPTI NLVRDLPQGFSALEPLVDLP IGINITRFQTLALHRSYL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

