

Human IL-12B protein (Recombinant) (C-His) (STJP001280)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3593
Gene Symbol	IL12B
UniProt ID	IL12B_HUMAN
Immunogen Region	Ile23-Ser328
Immunogen Sequence	IWELKKDVYVVELDWYPDAP GEMVVLTCDTPEEDGITWTL DQSSEVLGSGKTLTIQVKEF GDAGQYTCHKGGEVLSHSL LLHKKEDGIWSTDILKDQKE PKNKTLRCEAKNYSGRFTC WWLTTISTDLTFSVKSSRGS SDPQGVTCGAATLSAERVRG DNKEYEYSVECQEDSACPA EESLPIEVMVDVHKLKYEN YTSSFFIRDIKPDPPKNLQ LKPLKNSRQVEVSWEYPDT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

