

SARS-CoV-2 Spike RBD protein (Active Recombinant) (C-His&Avi) (STJP001282)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4 or Supplied as a 0.22 Mu m filtered solution in PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Arg319-Phe541
Immunogen Sequence	RVQPTESIVRFPNITNLCPF GEVFNATRFASVYAWNRRKRI SNCVADYSVLVNSASFSTFK CYGVSPTKLNDLCFTNVYAD SFVIRGDEV RQIAPGQTGKI ADYNYKLPDDFTGCVIAWNS NNLDSKVGNGYNYLYRFRK SNLKPFRDISTEIQAGST PCNGVEGFNCYFPLQSYGFQ PTNGVGYQPYRVVLSFELL HAPATVCGPKKSTNLVKNKC VNF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

