

Mouse Siglec-15/CD33L3 protein (Recombinant) (C-His) (STJP001283)

GENERAL INFORMATION

| | |
|---------------|--------------|
| Product Type | Proteins |
| Host / Source | HEK293 cells |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Formulation | Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4. |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|--------------------|---|
| Immunogen Region | Arg24-Thr262 |
| Immunogen Sequence | RRDASGDLLNTEAHSAPAQR WSMQVPAEVNAEAGDAAVLP CTFTHPHRHYDGPLTAIWRS GEPYAGPQVFRCTAAPGSEL CQTALSLHGRFRLLGNPRRN DLSLRVERLALADSGRYFCR VEFTGDAHDRYESRHRGVRLR VTAAAPRIVNISVLPGPAHA FRALCTAEGEPPPALAWSGP APGNSSAALQGQGHGYQVTA ELPALTRDGRYTCTAANSLG RAEASVYLFRRFHGAPGTST |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|

