

Mouse Siglec-15/CD33L3 protein (Recombinant) (C-hFc) (STJP001284)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Arg24-Thr262
Immunogen Sequence	RRDASGDLLNTEAHSAPAQR WSMQVPAEVNAEAGDAAVLP CTFTHPHRHYDGPLTAIWRS GEPYAGPQVFRCTAAPGSEL CQTALSLHGRFRLLGNPRRN DLSLRVERLALADSGRYFCR VEFTGDAHDRYESRHRGVRLR VTAAAPRIVNISVLPGPAHA FRALCTAEGEPPPALAWSGP APGNSSAALQGQGHGYQVTA ELPALTRDGRYTCTAANSLG RAEASVYLFRRFHGAPGTST

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

