

Cynomolgus Siglec-15/CD33L3 protein (Recombinant) (C-hFc&His) (STJP001286)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, 300mM NaCl, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Phe20-Thr263
Immunogen Sequence	FVRTKIDTTENLLNTEVHSS PAQRWSMQVPAEVSAAAGDA AVLPCFTFTHPHRHYDGPLTA IWRAGEPYAGPQVFRCAAR GSELCQTALSLHGRFRLGN PRRNDLSLRVERLALADDDR YFCRVEFAGDVHDDRYESRHG VRLHVTAAPRIINISVLPGP AHAFRALCTAEGEPPPALAW SGPALGNSSAAVPSSGQGHG HLVTAELPALNHDGRYTCTA ANSLGRSEASVYLFRFHGA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

