

SARS-CoV-2 Spike S1+S2 ECD (S-ECD) (D614G) protein (Active Recombinant) (C-His) (STJP001290)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Val11-Gln1208Asp614Gly
Immunogen Sequence	VSSQCVNLTTRTQLPPAYTN SFTRGVYYPDKVFRSSVLHS TQDLFLPFFSNVTWFHAIHV SGTNGTKRFDNPVLPFNDGV YFASTEKSNIRGWIFGTLL DSKTQSLIVNNAITNVVIKV CEFQFCNDPFLGVYYHKNNK SWMESEFRVYSSANNCTFEY VSQPFLMDLEGKQGNFKNLR EFVFKNIDGYFKIYSKHTPI NLVRDLPPQGFSALEPLVDLP IGINITRFQTLALHRSYL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

